

IV. NOTES ON USE

1. Use mild detergents for cleaning and disinfection. Do not use corrosive or abrasive products!
NOTE: Do not use agents containing naphtha, solvents, acids, chlorine etc. for cleaning soft surfaces!
2. Product marking is performed in the form of an identification label attached to each device. The label contains the following information:
 - number indicating the date of manufacture:
 - LOT: RRRR-MM (RRRR means the year, MM means the month, e.g.: 2023-08)
 - name and type of product e.g.: REF: KS/MR/Sk
 - allowed body weight of the user
3. **All devices (except for caps) have a 2-year warranty.**
4. **Use the devices as directed!** Any bending of the metal parts with visible cracks on the surface exclude the device from further use.
Check the device before and during use to ensure that it has no scratches, cracks, etc.
5. The production service during the warranty period (24 months) includes replacement of parts or products, which are found to have any manufacturing defects.
6. The MIKIRAD operational service functions directly at the points of sale of the products.
7. MIKIRAD guarantees the production of spare parts for 10 years from the date of manufacture of a product.

INSTRUCTIONS FOR USE COMMUNE CHAIRS

MIKIRAD manufactures **commode chairs** as medical products designed for persons with motor disabilities, classified as restroom aids according to the 09 12 class of the PN-EN ISO 9999 standard.

PN EN ISO 9999	NAME	TYPE
09 12 03	- Folded commode chair - Wheelless - Wheeled with 4 wheels (2 brakes)	KS/MR/Sk KS/MR/Sk.4k

I. BACKGROUND

Commode chairs facilitate the use of restrooms by persons with motor disabilities.

Our standard products are designed for users up to 120 kg in weight.




The chair is made from top-quality light and strong aluminium tubes joined together with plastic connectors.

Each chair has a removable toilet seat bolted to the rear tube of the structure and a container with a tight cover made from an easy-to-clean plastic. The legs of the chair are secured with non-slip soft plastic ferrules or provided with 4 wheels, including 2 brakes. The palm supports are padded with soft polyurethane foam.

Dimensions of a complete set:

- Seat depth: 43 cm
- Seat width: 48 cm
- Weight: 4,46/4,20 kg (with / without wheels)

The height of the chair can be adjusted within the 530-630 mm range in 20 mm increments by extending the adjusting tubes from the chair body and engaging the latches in the tube holes.

 Mikirad Sp. z o.o. 26-600 Radom, ul. Małopolska 16 tel.: +48 48 383 08 85 +48 48 383 08 87		 System zarządzania ISO 9001:2015 www.tuv.com ID 9105018398
---	---	---



Picture 1. Folded commode chair wheelless



Picture 2. Folded commode chair with wheels

II. OPERATION

1. Raise the toilet seat and slide in the container from the front and interlock it with the hangers placed on the supports under the seat.
To facilitate cleaning, pour some water to the container (just enough to cover its bottom) before using the chair.
2. After installing the container, remove the container's cover and lower the toilet seat.
3. To remove the container, raise the seat, close the container tightly with the cover and carefully pull the container to disengage it from the hangers.
4. Clean the container before re-use.

III. PROCEDURE FOR STEPWISE ADJUSTMENT OF PRODUCT HEIGHT (Fig. 1)

1. Undo (loosen) the nut of the damping sleeve.
2. Push the double latch out of the holes by pushing both fronts of the latch wings towards the back using two thumbs, so as to move the wings and slide the steel pegs out of the rod holes.
Do not bend the latch wings as they may break or lose elasticity.
3. Slide the adjustment rod along the axis to the appropriate height for the user; remember that the holes in the adjustment rod must ideally overlap the holes in the shaft.
4. Insert the steel pegs of the latch into the correctly positioned holes and push the latch using your hand until the inside surface of the latch is in full contact with the rod shaft.
5. Slightly tighten the nut of the damping sleeve.




Fig.1 Stepwise adjustment

IV. UNFOLDING AND FOLDING A COMMODE CHAIR

To unfold the chair, position it so that the manufacturer's label placed on the bent tube of the middle wing faces forward. Spread the two side wings and connect their articulated tubes so that the two tubes are locked once the thinner tube is inserted into the thicker one and the latch is engaged. Then install the rack supporting the container and the seat by pressing it into the parallel lateral tubes of the side wings. Slide the container into the wire hangers fitted near the rack. Reverse the foregoing procedure to fold the walker. After disengaging the latch, spread apart the wings to disconnect the tubes. Swing one of the wings forward and the other wing backward to flatten the components as much as possible.